


A 3D printed interbody material featuring a nano-rough surface showing deeper implant osseointegration with antibacterial characteristics uniquely combined in one radiolucent bone-like material.¹

Why compromise?

- Participates in fusion¹
- Bone-like mechanical properties¹
- Radiolucent
- Antibacterial characteristics¹



Fortilink[®]-TS

IBF System

with TETRAfuse[®]3D
TECHNOLOGY 

- 22, 26 and 32mm lengths
- Heights from 8-17mm (1mm increments)
- Lordotic angle: parallel and 6°
- Self distracting bulleted nose
- Patterned graft window designed to aid in bone-to-implant osseointegration



Sanderson's stained histology image (100x) from ovine study indicating bone cells (pink) interdigitating with TETRAfuse sample (tan)

Learn more: www.TETRAfuse3D.com

1. Data on file at RTI Surgical, Inc. Animal and In vitro data may not be representative of clinical experience.

Please refer to the labeling for clinical applications, warnings, precautions and other instructions for use.



RTI Surgical, Inc.

® indicates U.S. trademark registration. All trademarks and/or images are the property of their respective owners or holders.

©2017 RTI Surgical, Inc. All rights reserved.

Distributed by:

RTI Surgical, Inc.

375 River Park Circle
Marquette, MI 49855 USA

t: 800.557.9909

www.rtisurgical.com

10197 R0 11/2017