



nanOss Strata™ Advanced Bone Graft Powered by Vantage™ Technology



Two configurations. One advanced
technology for optimal bone healing.

25
YEARS

Building a Legacy
of Excellence


Two Configurations, One Technology

nanOss Strata™ products feature Vantage™ Technology — a biomimetic, nanocrystalline structure that supports cell attachment and reliable handling for diverse grafting needs.



nanOss Strata™ Moldable

- Moldable paste with collagen
- 1–2 mm Vantage™ granules
- Absorption: 1–2 days; resorption: 6–12 mo
- Firm, cohesive, resists irrigation



nanOss Strata™ Strip

- Collagen-enhanced strip
- Rapidly hydrates, flexible, easy to cut
- 1–2 mm Vantage™ granules

Product Code	Description	Size
450615	nanOss Strata moldable, 1.5 cc	1.5 cc
450603	nanOss Strata moldable, 3 cc	3 cc
450606	nanOss Strata moldable, 6 cc	6 cc
450612	nanOss Strata moldable, 12 cc	12 cc

Product Code	Description	Size
450508	nanOss Strata Strip, 50x25x8mm	50x25x8mm
450108	nanOss Strata Strip, 100x25x8mm	100x25x8mm
450506	nanOss Strata Strip, 50x25x6mm	50x25x6mm
450176	nanOss Strata Strip, 50x17x6mm	50x17x6mm



888.886.9354
CS@xtantmedical.com
xtantmedical.com